



DURAN® Gaswaschflaschen-Unterteile, nach Drechsel

DURAN® Gas washing bottles Lower Part, Drechsel type

	Inhalt [ml] Capacity [ml]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	100	10	24714240
	250	10	24714360
	500	10	24714440 ¹

- Mit mit NS 29/32, Bauhöhe 200 mm
- Nach DIN ISO 12463
- ¹ Ohne Fuß, ergänzend zur DIN
- Temperaturbeständigkeit: + 200°C - nicht isoliert, Minustemperatur limitiert durch Schliff fett/dichtung
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3
- Zusätzlich bei Druck: beachten Sie die Druckbeständigkeit von Klammern, Gabelklammern und Verschraubungen für NS-System

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- With NS 29/32, height 200 mm
- According to DIN ISO 12463
- ¹ Without base, supplementary to DIN
- Temperature resistance: + 200°C - not insulated, minus temperature limited by grinding grease/sealing
- Pressure resistance: -1/+0.5 barg
- Material in contact with product: Borosilicate glass 3.3
- Additionally for pressure: note the pressure resistance of clamps, fork clamps and screw connections for NS system

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.