



DURAN® Mikro-Filterkerzen mit engem Rohr

DURAN® Micro filter candles with narrow tube

	Porosität Porosity	Kerze Cartridge		d [mm]	h [mm]	Verpackungseinheit/ Stück Pack/Quantity	Artikelnr. Article no.
		d ₁ [mm]	h ₁ [mm]				
DWK	1-4	13	25	8	200	10	258573X ¹

- ¹ In die Best.-Nr. ist statt X die gewünschte Porosität einzusetzen, siehe auch Seite 240.
- Temperaturbeständigkeit: + 200°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- ¹ The required porosity should be inserted in place of X in the catalogue number, also see page 240.
- Temperature resistance: + 200°C
- Pressure resistance: -1/+0,5 barg
- Material in contact with product: Borosilicate glass 3.3

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.