



## DURAN® Schlitzsieb-Platten

## DURAN® Slit sieves

	Platten-Ø [mm] Disc. dia. [mm]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	48	10	21340310 <sup>1</sup>
	60	10	21340340
	73	10	21340390
	95	5	21340440
	120	1	21340510

- <sup>1</sup> Einsetzbar im DURAN®-Schraubfiltergerät, Best.-Nr. 24 720 50, als Träger für Membran- und Papierfilter.
- Temperaturbeständigkeit: + 200°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- <sup>1</sup> Can be used in DURAN® screwed filter apparatus, catalogue no. 24 720 50, as support for membrane and paper filters.
- Temperature resistance: + 200°C
- Pressure resistance: -1/+0,5 barg
- Material in contact with product: Borosilicate glass 3.3

*Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.*