



**Einhängeverbinder, Gestellbau,
LF Laborplanflansch**

*Support device, Scaffolding, LF
laboratory flat flange*

Für Nennweite DN <i>For nominal size DN</i>	Spannschrauben <i>Clamping screws</i>	Artikelnr. <i>Article no.</i> 13 mm Ø	Artikelnr. <i>Article no.</i> 21 mm Ø	Artikelnr. <i>Article no.</i> 27 mm Ø	Artikelnr. <i>Article no.</i> 30 mm Ø
60	4 x M5	CH060/13-L	CH060/21-L	CH060/27-L	-
100	4 x M6	CH100/13-L	CH100/21-L	CH100/27-L	CH100/30-L
120	4 x M6	CH120/13-L	CH120/21-L	CH120/27-L	CH120/30-L
150	4 x M6	CH150/13-L	CH150/21-L	CH150/27-L	CH150/30-L
200	4 x M6	CH200/13-L	CH200/21-L	CH200/27-L	CH200/30-L

- Massive Ausführung
- Komplett mit 2 Tragebügeln für Gestell 3/4" / 27 mm
- Achtung: bei Gefäßen kurzer Bauform werden breitere Tragebügel benötigt
- Material: Edelstahl 1.4301 (AISI 304)

- *Solid construction*
- *Complete with 2 carrying handles for rack 3/4" / 27 mm*
- *Attention: for vessels of short design wider support brackets are needed*
- *Material: stainless steel 1.4301 (AISI 304)*