



**GL-DURAN® Kunststoff-Oli-  
ven, gebogen, aus PP<sup>1</sup>**

**GL-DURAN® Plastic hose  
connections, bent, from PP<sup>1</sup>**

	Für Nenngröße For nominal size	AD [mm] OD [mm]	ID [mm] ID [mm]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	GL 14	8,6	4,0	10	29247050
	GL 18	10,8	6,8	10	28621260

- Material Olive: PP
- Material Dichtring: Silikon
- Temperaturbeständigkeit: -50/+150°C
- Druckbeständigkeit: -1/+3 barg
- Material produktberührend: PP

- Material olive: PP
- Material sealing ring: silicone
- Temperature resistance: -50/+150°C
- Pressure resistance: -1/+3 barg
- Material in contact with product: PP