



## GL-Dichtungen mit Bohrung

## GL seals with bore

	Nenngröße GL Nominal size GL	Größe d x d <sup>1</sup> [mm] Size d x d <sup>1</sup> [mm]	Für Rohr Ø [mm] For tube Ø [mm]	Bohrung D [mm] Bore D [mm]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
BOLA	14	12 x 6	5,5-6,5	6	10	CSPB01406
	18	16 x 6	5,5-6,5	6	10	CSPB01806
	18	16 x 8	7,5-9,0	8	10	CSPB01808
	18	16 x 10	9,0-11,0	10	10	CSPB01810
	25	22 x 8	7,5-9,0	8	10	CSPB02508
	25	22 x 10	9,0-11,0	10	10	CSPB02510
	25	22 x 12	11,0-13,0	12	10	CSPB02512
	32	29 x 10	9,0-11,0	10	10	CSPB03210
	32	29 x 12	11,0-13,0	12	10	CSPB03212
	32	29 x 14	13,0-15,0	14	10	CSPB03214
	32	29 x 16	15,0-17,0	16	10	CSPB03216
	32	29 x 18	17,0-19,0	18	10	CSPB03218
	45	42 x 26	25,0-27,0	26	10	CSPB04526
	45	42 x 32	31,0-33,0	32	10	CSPB04532

- Material: Silikon PTFE
- Temperaturbeständigkeit: -60/+230°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: PTFE

- Material: silicone/PTFE
- Temperature resistance: -60/+230°C
- Pressure resistance: -1/+0,5 barg
- Material in contact with product: PTFE