



GL-Schraubkappen, geschlossen

GL screw caps, closed

	Nenngröße GL Nominal size GL	Temperaturbeständigkeit [°C] temperature resistance [°C]	Verpackungseinheit/ Stück Pack/Quantity	Artikelnr. Article no.
DWK	14	140	10	CS014
	18	140	10	CS018
	25	140	10	CS025
	32	140	10	CS032
	45	140	10	CS045
BOLA	14	230	10	CS014HT
	18	230	10	CS018HT
	25	230	10	CS025HT
	32	230	10	CS032HT
	45	230	10	CS045HT

- Material 140°C: PBT
- Temperaturbeständigkeit: -50/+140°C
- Material 230°C: PPS
- Temperaturbeständigkeit: -20/+230°C, HT-Variante, drucklos, max 200°C bei 0,5 barg
- Material Flachdichtung: Silikon/PTFE
- Druckbeständigkeit: -1/+10 barg
- Material produktberührend: PTFE

- *Material 140°C: PBT*
- *Temperature resistance: -50/+140°C*
- *Material 230°C: PPS*
- *Temperature resistance: -20/+230°C. HT version, pressureless, max 200°C at 0,5 barg*
- *Material flat gasket: silicone/PTFE*
- *Pressure resistance: -1/+10 barg*
- *Material in contact with product: PTFE*