



GL-Schraubkappen, offen

GL screw caps, open

	Nenngröße GL Nominal size GL	Bohrung D [mm] Bore D [mm]	Temperaturbeständigkeit [°C] Temperature resistance [°C]	Stk./Verpackungseinheit Piece/packing unit	Artikelnr. Article no.
NORMAG	14	9,5	140	10	CSO014
	18	11	140	10	CSO018
	25	15	140	10	CSO025
	32	20	140	10	CSO032
	45	34	140	10	CSO045
	14	9,5	250	10	CSO014HT
	18	11	250	10	CSO018HT
	25	15	250	10	CSO025HT
	32	20	250	10	CSO032HT
	45	34	250	10	CSO045HT

- Material 140°C: PBT
- Temperaturbeständigkeit: -50/+140°C
- Material 250°C: PPS
- Temperaturbeständigkeit: -20/+250°C. HT-Variante, drucklos, max 200°C bei 0,5 barg
- Druckbeständigkeit: -1/+10 barg
- Nicht produktberührend

- Material 140°C: PBT
- Temperature resistance: -50/+140°C
- Material 250°C: PPS
- Temperature resistance: -20/+250°C. HT variant, unpressurized, max 200°C at 0.5 barg
- Pressure resistance: -1/+10 barg
- Not in contact with product