


NORMAG



HF-Außenzentrierringe mit FKM Dichtung

HF External centering rings with FKM seal

	Nenngröße HF Nominal size HF	D [mm] D [mm]	H [mm]	Artikelnr. Article no.
NORMAG	10/16	32	5	CSVHF16
	25	42	5	CSVHF25
	40	57	5	CSVHF40
	50	77	5	CSVHF50

- Material: Aluminium/FKM
- Temperaturbeständigkeit: -40/+180°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: FKM

- *Material: aluminium/FKM*
- *Temperature resistance: -40/+180°C*
- *Pressure resistance: -1/+0.5 barg*
- *Material in contact with product: FKM*