


NORMAG



HF-Flanschverbindungen

HF Flange connections

	Nenngröße HF Nominal size HF	D [mm] D [mm]	Artikelnr. Article no.
NORMAG	10/16	55	CAHF16
	25	68	CAHF25
	40	85	CAHF40
	50	105	CAHF50

- Material: Aluminium eloxiert natur
- Temperaturbeständigkeit: -60/+300°C
- Druckbeständigkeit: -1/+0,5 barg
- Nicht produktberührend

- *Material: anodised aluminium natural*
- *Temperature resistance: -60/+300°C*
- *Pressure resistance: -1/+0.5 barg*
- *Non-product-contacting*