


**NORMAG**



**KF/PF-Übergangsdichtung,  
universell**

*KF/PF transition seal, universal*

	Nennweite DN 1/2 Nominal size DN 1/2	L [mm]	Artikelnr. Article no.
NORMAG	15	6	CGE015
	25	7	CGE025
	40	8	CGE040
	50	8	CGE050

- Material: PTFE
- Temperaturbeständigkeit: -100/+200°C
- Druckbeständigkeit: -1/ max. +6 barg (abhängig von DN)
- Material produktberührend: PTFE

- *Material: PTFE*
- *Temperature resistance: -100/+200°C*
- *Pressure resistance: -1/ max. +6 barg (depending on DN)*
- *Material in contact with product: PTFE*