


NORMAG



KF-PTFE Universaldichtungen

KF universal seals PTFE

	Nennweite DN Nominal size DN	Artikelnr. Article no.
NORMAG	15	CGR015-K
	25	CGR025-K
	40	CGR040-K
	50	CGR050-K
	80	CGR080-K
	100	CGR100-K
	150	CGR150-K
	200	CGR200
	300	CGR300

- Material: PTFE
- Temperaturbeständigkeit: -100/+200°C
- Druckbeständigkeit: -1/ max. +6 barg (abhängig von DN)
- Material produktberührend: PTFE

- *Material: PTFE*
- *Temperature resistance: -100/+200°C*
- *Pressure resistance: -1/ max. +6 barg (depending on DN)*
- *Material in contact with product: PTFE*