

**NORMAG**



**LF-Dichtungen Silikon,  
FEP ummantelt**

*LF seals silicone, FEP coated*

	Nenngröße DN <i>Nominal size DN</i>	Rohr D x s [mm] <i>Tube D x s [mm]</i>	Artikelnr. <i>Article no.</i>
NORMAG	60	75 x 4	CSFLF060
	100	110 x 4	CSFLF100
	120	133 x 4	CSFLF120
	150	157 x 5	CSFLF150
	200	215 x 5	CSFLF200

- Material: Silikon/FEP ummantelt
- Temperaturbeständigkeit: -100/+200°C
- Druckbeständigkeit: -1/+2,5 barg
- Material produktberührend: FEP

- *Material: Silicone/FEP coated*
- *Temperature resistance: -100/+200°C*
- *Pressure resistance: -1/+2.5 barg*
- *Material in contact with product: FEP*