

**NORMAG**



**NS-Gabelklammern**

*NS fork clamps*

	Nenngröße NS Nominal size NS	Artikelnr. Article no.
NORMAG	14	CFNS14
	19	CFNS19
	29	CFNS29
	45	CFNS45

- Material: Edelstahl
- Temperaturbeständigkeit: -60/+300°C
- Druckbeständigkeit: -1/+0,1 barg
- Nicht produktberührend

- *Material: stainless steel*
- *Temperature resistance: -60/+300°C*
- *Pressure resistance: -1/+0.1 barg*
- *Non-contact with product*