


NORMAG



**PF-Schellenringverbindungen,
Glas-Glas**

*PF clamp ring connections,
glass-glass*

	Nennweite DN Nominal size DN	D [mm]	H [mm]	TK [mm]	Schrauben n x M x L [mm] Screws n x M x L [mm]	Artikelnr. Article no.
NORMAG	15	70	39	50	3 x M6 x 65	CP015-P
	25	90	61	70	3 x M8 x 90	CP025-P
	40	109	66	86	3 x M8 x 100	CP040-P
	50	122	73	98	3 x M8 x 105	CP050-P
	80	160	87	133	6 x M8 x 120	CP080-P
	100	204	98	178	6 x M8 x 130	CP100-P
	150	280	100	254	6 x 118 x 150	CP150-P

- Material Schellenringe: Duroplastisches Polyester mit Glas/Mineralanteil
- Material Einlage: Phenol-Formaldehyd
- Material Schrauben/Muttern/Scheiben: V4A
- Material Federn: Edelstahl 1.4310
- Temperaturbeständigkeit: -50/+200°C, nicht isoliert
- Druckbeständigkeit: -0,5/ max. +6 barg (abhängig von DN)
- Nicht produktberührend

- *Material clamp rings: Duroplastic polyester with glass/mineral content*
- *Material insert: Phenol-Formaldehyde*
- *Material screws/nuts/washers: V4A*
- *Material springs: stainless steel 1.4310*
- *Temperature resistance: -50/+200°C, not insulated*
- *Pressure resistance: -1/ max. +6 barg (depending on DN)*
- *Non-product-contacting*