

NORMAG



S/KS/RS-Gabelklammern

S/KS/RS-fork clamps

	Nenngröße S Nominal size S	Nenngröße KS Nominal size KS	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
NORMAG	13	12	5	CFS13
	19	18	5	CFS19
	29	28	5	CFS29
	35	35	2	CFS35
	40	-	2	CFS40
	51	50	1	CFS51
	64	-	1	CFS64

- Material: Edelstahl
- Temperaturbeständigkeit: -60/+300°C
- Druckbeständigkeit: -1/+0,1 barg
- Nicht produktberührend

- *Material: stainless steel*
- *Temperature resistance: -60/+300°C*
- *Pressure resistance: -1/+0.1 barg*
- *Non-contact with product*