


NORMAG



S/KS/RS-Schellenring- verbindungen

S/KS/RS clamp ring connections

Nenngröße S/KS/RS Nominal size S/KS/RS		Artikelnr. Article no.
NORMAG	35	CCS35
	40	CCS40
	51	CCS51
	64	CCS64
	75	CCS75
	90	CCS90

Nenngröße S/KS/RS Nominal size S/KS/RS		Vakuumentel einseitig vacuum jacket one-sided	Artikelnr. Article no.
NORMAG	35		CCSJ35
	40		CCSJ40
	51		CCSJ51
	64		CCSJ64
	75		CCSJ75
	90		CCSJ90

Nenngröße S/KS/RS Nominal size S/KS/RS		Vakuumentel zweiseitig vacuum jacket two-sided	Artikelnr. Article no.
NORMAG	35		CCSJS35
	40		CCSJS40
	51		CCSJS51
	64		CCSJS64
	75		CCSJS75
	90		CCSJS90

- Material Schellenringe: Aluminium pulverbeschichtet
- Material Einlage: PTFE
- Material Schrauben/Mutter: V4A
- Material Feder: Edelstahl 1.4310
- Temperaturbeständigkeit: -60/+300°C
- Druckbeständigkeit: -1/max. +6 barg, in Abhängigkeit der Größe

- Material clamp rings: aluminium powder-coated
- Material insert: PTFE
- Material screws/nut: V4A
- Material spring: stainless steel 1.4310
- Temperature resistance: -60/+300°C
- Pressure resistance: -1/max. +6 barg, depending on the size