

NORMAG



SP-Dichtungen, FFKM

SP seals, FFKM

	Nenngröße SP Nominal size SP	Artikelnr. Article no.
NORMAG	13/5	CSKSP13
	19/9	CSKSP19
	29/15	CSKSP29
	35/20	CSKSP35
	41/25	CSKSP41
	64/40	CSKSP64
	75	CSKSP75
	90	CSKSP90

- Material: FFKM
- Temperaturbeständigkeit: -15/+270°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: FFKM

- *Material: FFKM*
- *Temperature resistance: -15/+270°C*
- *Pressure resistance: -1/+0.5 barg*
- *Material in contact with product: FFKM*