

# NORMAG



## SP-Dichtungen FKM

## SP seals FKM

	Nenngröße SP Nominal size SP	Stk./Verpackungseinheit Piece/packing unit	Artikelnr. Article no.
NORMAG	13/5	5	CSVSP13
	19/9	5	CSVSP19
	29/15	5	CSVSP29
	35/20	5	CSVSP35
	41/25	5	CSVSP41
	64/40	5	CSVSP64
	75	1	CSVSP75
	90	1	CSVSP90

- Material: FKM
- Temperaturbeständigkeit: -40/+180°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: FKM

- *Material: FKM*
- *Temperature resistance: -40/+180°C*
- *Pressure resistance: -1/+0.5 barg*
- *Material in contact with product: FKM*