


# NORMAG



## SP-Gabelklammern

## SP fork clamps

Verbindungselemente  
Joints and connectors

	Nenngröße SP Nominal size SP	Artikelnr. Article no.
NORMAG	13	CFS13
	19	CFS19
	29	CFS29
	35	CFS35
	40	CFS40
	64	CFS64

- Material: Edelstahl
- Temperaturbeständigkeit: -60/+300°C
- Druckbeständigkeit: -1/+0,1 barg
- Nicht produktberührend

- *Material: stainless steel*
- *Temperature resistance: -60/+300°C*
- *Pressure resistance: -1/+0.1 barg*
- *Non-contact with product*