



Security Flex Dichtungen

Security flex seals

Rohr AD [mm] Tube AD [mm]	Gewinde GL 14 Artikelnr. Thread GL 14 Article no.	Gewinde GL 18 Artikelnr. Thread GL 18 Article no.	Gewinde GL 25 Artikelnr. Thread GL 25 Article no.	Gewinde GL 32 Artikelnr. Thread GL 32 Article no.	Gewinde GL 45 Artikelnr. Thread GL 45 Article no.
0	CSF014000 ¹	CSF018000 ¹	CSF025000 ¹	CSF032000 ¹	CSF045000 ¹
1,5	CSF014015	CSF018015	CSF025015	-	-
1,6 (1/16")	CSF014016	CSF018016	-	-	-
1,9	CSF014019	CSF018019	-	-	-
2	CSF014020	CSF018020	CSF025020	-	-
3	CSF014030	CSF018030	CSF025030	-	-
3,2 (1/8")	CSF014032	CSF018032	CSF025032	-	-
4	CSF014040	CSF018040	CSF025040	CSF032040	-
6	CSF014060	CSF018060	CSF025060	CSF032060	CSF045060
6,4 (1/4")	CSF014064	CSF018064	CSF025064	CSF032064	-
8	CSF014080	CSF018080	CSF025080	CSF032080	CSF045080
9,5	-	CSF018095	CSF025095	CSF032095	-
10	CSF014100	CSF018100	CSF025100	CSF032100	CSF045100
12	-	-	CSF025120	CSF032120	CSF045120
14	-	-	CSF025140	CSF032140	CSF045140
16	-	-	-	CSF032160	CSF045160
18	-	-	-	CSF032180	CSF045180
20	-	-	-	CSF032200	CSF045200
22	-	-	-	-	CSF045220
25,4 (1")	-	-	-	-	CSF045254
26	-	-	-	-	CSF045260
30	-	-	-	-	CSF045300
32	-	-	-	-	CSF045320

NORMAG

NORMAG

- Flexible hoch abdichtende Verbindung, rutschfest
 - Verhindert selbstständiges lösen äußerst sicher
 - Auch bei größeren Temperaturschwankungen dicht
 - ¹Blindkappe ohne Bohrung
 - Grüne Dichtung für Zoll-Maße
 - Transparente Dichtung für mm-Maße
 - Für Vakuumeinsatz geeignet
 - Material: Silikon/PTFE
 - Spezielles Silikon, dass Temperaturbelastungen lange stand hält, ohne weich oder klebrig zu werden
 - Temperaturbeständigkeit: -50/+250°C
 - Druckbeständigkeit: -1/+0,5 barg
 - Material produktberührend: PTFE
 - Sondermaße auf Anfrage
- *Flexible high sealing connection, non-slip*
 - *Prevents independent loosening extremely securely*
 - *Tight even with large temperature fluctuations*
 - *¹ blind cap without hole*
 - *Green seal for inch dimensions*
 - *Transparent seal for mm dimensions*
 - *Suitable for vacuum use*
 - *Material: silicone/PTFE*
 - *Special silicone that withstands temperature stresses for a long time without becoming soft or sticky*
 - *Temperature resistance: -50/+250°C*
 - *Pressure resistance: -1/+0,5 barg*
 - *Material in contact with product: PTFE*
 - *Special sizes on request*