



**Security Flex Multi Plus**

*Security flex multi plus*

Gewinde GL Thread GL		Artikelnr. Article no.
45		CMPP045

NORMAG

Ersatzteilset Spare part kit		
Bezeichnung Description	Artikelnr. Article no.	
Kappe geschlossen, schwarz/rot (Bild 1)	CMPPC	
Kappe offen, schwarz (Bild 1)	CMPCO	
PTFE Trägersystem, weiße Dichtung und GL045 Präzisionskappe (Bild 2)	CMPP045	
Security Flex Plus Dichtung	Seite 90	

NORMAG

- Flexible Möglichkeit zur sicheren Verbindung mehrerer Rohre, Schläuche, Sensoren, ...mit einem GL45
- Zur Abdichtung werden Security Flex Plus Dichtungen verwendet
- Nicht benötigte Anschlüsse können mit Blindstopfen versehen werden
- CMPP045 enthält die GL045 Präzisionskappe und das PTFE-Trägersystem, dazu werden die Kappen und Dichtungen noch benötigt
- Material: glasfaserverstärktes PPS
- Temperaturbeständigkeit: -20/+230°C
- Material produktberührend: -
- Druckbeständigkeit:-1/+6 barg
- Flexible way to securely connect multiple pipes, hoses, sensors, ...to a GL45.
- Security Flex Plus gaskets are used for sealing
- Unused connections can be fitted with blind caps
- CMPP045 includes GL045 precision cap and PTFE carrier system, caps and gaskets still required
- Material: glass fibre reinforced PPS
- Temperature resistance: -20/+230°C
- Material in contact with product: -
- Pressure resistance:-1/+6 barg