



**Security Flex Plus Dichtungen**    *Security flex plus seals*

Rohr AD [mm] Tube AD [mm]	Gewinde GL 14 Artikelnr. Thread GL 14 Article no.	Gewinde GL 18 Artikelnr. Thread GL 18 Article no.	Gewinde GL 25 Artikelnr. Thread GL 25 Article no.	Gewinde GL 32 Artikelnr. Thread GL 32 Article no.	Gewinde GL 45 Artikelnr. Thread GL 45 Article no.
0	CSFP014000 <sup>1</sup>	CSFP018000 <sup>1</sup>	CSFP025000 <sup>1</sup>	CSFP032000 <sup>1</sup>	CSFP045000 <sup>1</sup>
1,5	CSFP014015	CSFP018015	CSFP025015	-	-
1,6 (1/16")	CSFP014016	CSFP018016	-	-	-
1,9	CSFP014019	CSFP018019	-	-	-
2	CSFP014020	CSFP018020	CSFP025020	-	-
3	CSFP014030	CSFP018030	CSFP025030	-	-
3,2 (1/8")	CSFP014032	CSFP018032	CSFP025032	-	-
4	CSFP014040	CSFP018040	CSFP025040	CSFP032040	-
6	CSFP014060	CSFP018060	CSFP025060	CSFP032060	CSFP045060
6,4 (1/4")	CSFP014064	CSFP018064	CSFP025064	CSFP032064	-
8	CSFP014080	CSFP018080	CSFP025080	CSFP032080	CSFP045080
9,5	-	CSFP018095	CSFP025095	CSFP032095	-
10	CSFP014100	CSFP018100	CSFP025100	CSFP032100	CSFP045100
12	-	-	CSFP025120	CSFP032120	CSFP045120
14	-	-	CSFP025140	CSFP032140	CSFP045140
16	-	-	-	CSFP032160	CSFP045160
18	-	-	-	CSFP032180	CSFP045180
20	-	-	-	CSFP032200	CSFP045200
22	-	-	-	-	CSFP045220
25,4 (1")	-	-	-	-	CSFP045254
26	-	-	-	-	CSFP045260
30	-	-	-	-	CSFP045300
32	-	-	-	-	CSFP045320

NORMAG

# NORMAG

- Flexible hoch abdichtende Verbindung, rutschfest
  - Verhindert selbstständiges lösen äußerst sicher
  - Auch bei größeren Temperaturschwankungen dicht
  - Dichtung mit optimierter Haltefunktion für Schläuche und Sensoren
  - <sup>1</sup>Blindkappe ohne Bohrung
  - Maximale Dichteigenschaften für Schlauch/Rohrleitungen und GL-Gewinde
  - Transparente Dichtung für mm-Maße
  - Grüne Dichtung für Zoll-Maße
  - Material: Silikon/PTFE
  - Spezielles Silikon, dass Temperaturbelastungen lange stand hält, ohne weich oder klebrig zu werden
  - Temperaturbeständigkeit: -50/+250°C
  - Druckbeständigkeit: -1/+6 barg
  - Material produktberührend: PTFE
- *Flexible high sealing connection, non-slip*
  - *Prevents independent loosening extremely securely*
  - *Tight even with large temperature fluctuations*
  - *Seal with optimised holding function for hoses and sensors*
  - *<sup>1</sup>Blind cap without hole*
  - *Maximum sealing properties for hoses/pipes and GL threads*
  - *Transparent seal for mm-dimensions*
  - *Green seal for inch dimensions*
  - *Material: silicone/PTFE*
  - *Special silicone that withstands temperature stresses for a long time without becoming soft or sticky*
  - *Temperature resistance: -50/+250°C*
  - *Pressure resistance: -1/+6 barg*
  - *Material in contact with product: PTFE*