

**NORMAG**



**Security Plus Präzisionskappen**

*Security plus precision caps*

	Gewinde GL Thread GL	Artikelnr. Article no.
	18/14	CSPP018014 <sup>1</sup>
	14	CSPP014
	18	CSPP018
	25	CSPP025
	32	CSPP032
	45	CSPP045
NORMAG	18/14	CSPPEX018014 <sup>1,2</sup>
	14	CSPPEX014 <sup>2</sup>
	18	CSPPEX018 <sup>2</sup>
	25	CSPPEX025 <sup>2</sup>
	32	CSPPEX032 <sup>2</sup>
	45	CSPPEX045 <sup>2</sup>

- <sup>1</sup> in eine GL 18-Kappe eingeschraubtes GL 14-Gewinde, durch das ein 10mm-Rohr verbunden werden kann. Kappe kann weiterhin für GL18 verwendet werden
- <sup>2</sup> Ex-Ausführung
- Kappen aus Vollmaterial gedreht
- Hohe mechanische Belastbarkeit
- Material: glasfaserverstärktes PPS
- Temperaturbeständigkeit: -20/+230°C
- Material produktberührend: -
- Druckbeständigkeit:-1/+6 barg

- <sup>1</sup> GL 14 thread screwed into a GL 18 cap through which a 10mm tube can be connected. Cap can still be used for GL18
- <sup>2</sup> Explosion-proof design
- Caps turned from solid material
- High mechanical strength
- Material: glass fibre reinforced PPS
- Temperature resistance: -20/+230°C
- Material in contact with product: -
- Pressure resistance:-1/+6 barg