



Security Präzisionskappen

Security precision caps

	Gewinde GL Thread GL	Artikelnr. Article no.
NORMAG	18/14	CSP018014 ¹
	14	CSP014
	18	CSP018
	25	CSP025
	32	CSP032
	45	CSP045
	18/14	CSPEX018014 ^{1,2}
	14	CSPEX014 ²
	18	CSPEX018 ²
	25	CSPEX025 ²
	32	CSPEX032 ²
	45	CSPEX045 ²

	Bezeichnung Description	Artikelnr. Article no.
NORMAG	Adapter GL14 auf GL18 (beige) Adapter GL14 to GL18 (beige)	VSV018014
	Adapter GL14 auf GL18, leitfähig (schwarz) Adapter GL14 to GL18, conductive (black)	VSVEX018014

- ¹ in eine GL 18-Kappe eingeschraubtes GL 14-Gewinde, durch das ein 10mm-Rohr verbunden werden kann. Kappe kann weiterhin für GL18 verwendet werden
- ² Ex-Ausführung
- Kappen aus Vollmaterial gedreht
- Hohe mechanische Belastbarkeit
- Material: glasfaserverstärktes PPS
- Temperaturbeständigkeit: -20/+230°C
- Material produktberührend: -
- Druckbeständigkeit: -1/+0,5 barg

- ¹ GL 14 thread screwed into a GL 18 cap through which a 10mm tube can be connected. Cap can still be used for GL18
- ² Explosion-proof design
- Caps turned from solid material
- High mechanical strength
- Material: glass fibre reinforced PPS
- Temperature resistance: -20/+230°C
- Material in contact with product: -
- Pressure resistance: -1/+0.5 barg